






Quantitative data analysis




3/4/2022

Why?

- Example

Image quantification software


- IN Carta Image Analysis Software (Molecular devices)
- Nikon Image Analysis Software
- Image-Pro Plus
- Free Image Analysis Software



Fiji (ImageJ2)

ImageJ Docs

- Download
- Learn
- Extend
- Contribute
- Discuss
- Explore
- News
- Events
- Libraries
- Software
- NIH Image
- ImageJ
- ImageJ2
- Fiji
- More...



Fiji

Fiji is an image processing package—a “batteries-included” distribution of [ImageJ2](#), bundling a lot of plugins which facilitate scientific image analysis.

- **For users** - Fiji is [easy to install](#) and has an automatic update function, bundles a [lot of plugins](#) and offers comprehensive [documentation](#).
- **For developers** - Fiji is an open source project hosted in a [Git version control repository](#), with access to the source code of all internals, libraries and plugins, and eases the [development](#) and [scripting](#) of plugins.

Downloads

~ Download Fiji for your OS ~



64-bit



macOS



64-bit



32-bit



Other downloads



See the [Fiji downloads page](#) for Life-Line versions, etc.

License


Fiji is released as [open source](#) under the [GNU General Public License](#).

Fiji builds on top of the [ImageJ2](#) core, which is licensed under the permissive [BSD 2-Clause license](#).

Fiji also includes the original [ImageJ](#), which is [free of copyright](#).


Plugins and other components have [their own licenses](#).

Task 1. Basic measurement and Scale set




1. Set Scale
2. Measure Diameter, Area, and Distance between Nanoparticles (n=10)

Task 2. Fluorescent intensity using plot profile



1. Draw a Line within yellow boxed tissue region (gray matter)
2. Measure the Fluorescent intensity using plot profile ($n=10$ for each condition)
3. Quantitatively compare the intensity of each marker (NF160 for neurons and GFAP for astrocytes) between before and after injury
4. Statistically compare the results using P value (T-test): Significant difference between before and after injury?

Task 3. Protein intensity (Western blot)




Quantitatively compare the protein intensity (GAPDH) of each cell/tissue.

1. Select “Area of Interest”
2. Measure “Integrated density”


*Integrated density = Sum of Values of Pixels (similar to Area Under Curve)

Task 4. Counting Cells using a “Analyze Particle”





1. Drop “Task4 DAPI stained 1 8bit gray scale.tif”
2. Select the cells using a Threshold (click “dark background”)
3. Make it as a Binary image (under Process; Binary)
4. Apply “Watershed” for cell segmentation (under Process; Binary)
5. Count cells using a “Analyze Particle (click “outlines”)

Task 5. Nanoparticle size distribution analysis



1. Drop “Task5 Nanoparticle size distribution.jpg”
2. Set Scale
3. Select the nanoparticles using a Threshold (click “dark background”)
4. Apply “Watershed” for segmentation (under Process; Binary)
5. Count the nanoparticles using a “Analyze Particle (click “outlines”
6. Save the results in Excel for size distribution analysis

Task 6. Neurite tracking




Quantitatively compare the neurite outgrowth between Cortical neurons 1 and Cortical neurons 2

1. Set Scale
2. Select the neurites using a Threshold (click “dark background”)
3. Measure

Task 7. Stitching


13x18 Stitched Image (234 FOV)

Stitched image shows an entire well of a 6-well plate




Live cells labeled with Calcein, AM


Single FOV resolution




Tile 4




Tile 3




Tile 2



Tile 1



Stitched image



Stitch four tiles (Tiles 1-4) using a Grid/Collection Stitching (under Plugins-Stitching)

File names for tiles: **Task 7 Tile {ii}.tif**